


AJ18



C€0197 ISO13485:2003








www.ajaxdent.com

Patient Chair



- Danish motor from LINAK
 3 programmable position
 Reset position : easy to get on the chair

- Reset position: easy to get on the chair
 Last memory chair position: back to the
 previous treatment position easily
 FrgoMotion: combined movement of
 the backrest and the seat, prevents
 Stretching and compression of the
 patient's spine
 Chair-light concurrent on/off system
- Chair-ingth concurrent on/or system
 Chair-instrument lock: prevents chair moving when operating
 Seamless upholstery: easy to clean
 One master switch to control power, air and water



Chair backrest safety


Backrest and chair will be auto rised when they touch people's leg or foot hardly, ensure dentist's safety.





Articulate headrest

Any comfortable position for the patient's head can be fixed

Assistant element







AJ18 Configuration form

			AJ18V*Venturi)	AJ18D (Dry)	AJ18W (Wet)
			Hanging/Top	Hanging/Top	Hanging/Top
Dental light		FARO EDI	*	*	*
		FARO EDI with sensor	0	0	0
		ISA	0	0	0
		ISA with sensor	0	0	0
		Light II	0	2	-
		Light II with sensor	0	-	
Patient chair		Denish motor LINAK	*	*	
		Chinese motor	0		
		Articulated headrest	*	*	*
		Armrest	0	0	0
		3 programable position	*	*	*
		Reset postion and Last memory position			
		ErgoMotion			
		Chair safety system	*	*	*
		Chair-instrument lock		*	
Water unit		Moveable glass spittoon	1		
		Programmable cup filler and bowl rinsing	1 :		
		Suction filter		*	*
		Water heater		*	*
		Bottle for pure water	*	•	
		Air venturi system			
		Spittoon valve (Metasys)	-	-	
		Place selection valve (Metasys)	-	-	
		Amalgam separator (Metasys compact)	4	*	101
Air pressure hygiene system		Bottle for disinfectant	*	*	*
		Sanitation hood	*	*	
Delivery system	Position1	Three way syringe	*	*	
	Position2	Air tubine hose	О	-	-
		Air tubine hose with light	*	*	*
	position3	Air tubine hose	0	-	
		Air tubine hose with light	*	-	-
		NSK eletric motor with light	О	*	
		Bien air eletric motor MC2/MC3 with light	0	0	0
	position4	Air motor hose	0		
	p a a state of the	Air motor hose with light	*		
		NSK eletric motor with light	О	*	
		Bien air eletric motor MC2/MC3 with light	0	o	0
	position5	Woodpecker scaler	0	-	-
	positions	EMS scaler	*		
		Main control panel with 3 memories			
		Build-in X-ray viewer			*
Assitant element	Position1	Three way syringe		*	
	Position2	HVE with adjustable switch		*	*
	Position2 Position3		*	*	
		Saliva ejector with adjustable switch			
	Position4	Curing light		:	
		Assitant control panel			
		German suction hose	*	*	*
Foot control		Multi-function foot control			
Dental stool		Dental stool II (bigger)	*	*	*
		Another dental stool II	0	О	0
Intra oral camera		WS 619	0	-	-
		MiniCAM from Taiwan	*	*	*
		17" monitor		*	*

* Standard configuration o Optional - N/A

